


Návod k použití a technický popis
User manual and technical description
Gebrauchsanweisung und technische Beschreibung
Manual de usuario y descripción técnica


EffectaCare 10

EffectaCare 20


PrimaCare 10


PrimaCare 20

CliniCare 10

CliniCare 10

CliniCare 30

Pasivní matrace
Passive Mattresses
Passive Matratzen
Colchones estáticos

Kontakt / Contact / Kontakt / Contacto:

Linét spol. s r. o.
Želevčice 5,
274 01 Slaný
Czech Republic

Tel.: +420 312 576 111
Fax.: +420 312 522 668
E-Mail: info@linet.cz
Web: www.linet.com

Pasivní matrace / Passive mattresses / Passive Matratzen / Colchones estáticos

EffectaCare 10, EffectaCare 20, PrimaCare 10, PrimaCare 20, CliniCare 10, CliniCare 20, CliniCare 30

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All trademarks and brands are the property of the appropriate owners. The manufacturer reserves the right to changes in the contents of this manual that relate to the product's technical regulations. It is for this reason that the contents of this manual may indicate differences from the current manufacture of the product.

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ČEŠTINA

User manual and technical description

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1.1 Warning Notices

1.1.1 Types of Warning Notices

Warning notices are differentiated by the type of danger using the following key words:

- Caution - risk of material damage.
- Warning - risk of physical injury.
- Danger - risk of fatal injury.

1.1.2 Structure of Warning Notices

⚠ SIGNAL WORDS!	
Type and source of danger!	
⇒	Measures to avoid the danger.

1.2 Other Symbols

1.2.1 Instructions

Structure of instructions:

- ❖ Perform this step.
Result, if necessary.



1.2.2 Lists

Structure of bulleted lists:

- List level 1
 - List level 2

2 Used Symbols and Labels

2.1 Symbols and Labels on the Product

	Read user manual
	Attention
	CE mark of conformity with EU regulation
	Machine wash at 71 deg. Centigrade for 3 minutes
	Do not bleach.
	Tumble dry on low heat setting.
	Do not wring.
	Do not iron.
	Gentle dry cleaning with PCE (perchloroethylen).
	Do not use a cleaner containing Phenol.
	Cover Materials are Fire Retardent to BS7175, Sources 0, 1 and 5
	Regularly inspect the inside of the cover for contamination.
	Handwash with detergent. Initial temperature of hot water should not exceed 50 deg. Centigrade.
	Disinfect using solution containing $\leq 1,000\text{ppm}$ of Chlorine (see Cleaning/Disinfection).
	Rinse with water.
	Dry.
	Weight limit of mattress.
	Symbol for turning and rotating the mattress. See chapter "Using the Mattress."
	Identification and wash label.

3 Safety and Dangers

3.1 Safety Instructions

Use:

- ❖ It is necessary to read the user manual before using the mattress and perform maintenance according to the user manual.
- ❖ Follow the instructions carefully.
- ❖ Do not use a damaged mattress.
- ❖ If possible, check the mattress daily or between changing patients.
- ❖ Before positioning the patient the user must verify the condition of the mattress.
- ❖ Single sided mattresses are only to be used with the screen printing uppermost, and, if necessary, to check the mattress core has the profiling uppermost.
- ❖ The mattress must not be used if failures have been found that can do harm to the patient or mattress.
- ❖ The selection of the correct mattress type and application of the mattress must always be consulted with qualified clinical professional, especially in the case of an injury of the head, spinal cord or backbone.
- ❖ A patient with the indication of a higher risk of pressure ulcers than the mattress is designed for must not be placed on the mattress.
- ❖ Anti-decubitus mattresses are designed for the prevention and treatment of decubitus and this is why it is necessary to monitor the condition of the patient and his/hers skin continuously and to change his/hers position to prevent stressing risk and affected parts of his/hers body.
- ❖ The mattress must not be laid on the mattress platform if there is a danger the mattress could be disturbed or the cover torn with sharp parts of the bed or other sharp objects.
- ❖ During using of the mattress it is necessary to monitor the condition of the mattress core to ensure that it is not permanently compressed which will reduce the anti-decubitus properties of the product.
- ❖ Prolonged exposure of light can change the color tone of the mattress core. This change has no effect on physical properties of the mattress core. Because the polyurethane foams can have different colours. This has no effect on physical properties of the mattress core.
- ❖ The mattress cover may shrink slightly or its colour may fade. The mattress is still suitable for use provided the cover does not become a very tight fit (i.e can still bend under gentle pressure) and shows no signs of damage to the top coating.
- ❖ All mattresses with PU covers are fully waterproof and must not be covered with an additional waterproof sheet as this will stop its moisture vapour permeability function of mattress cover.
- ❖ Do not put the bed into Trendelenburg or Anti-Trendelenburg position, when the bed-ends are removed unless the bed frame is fitted with a retaining method to prevent the mattress from sliding up/down the bed when it at such angles.
- ❖ The mattress must be placed on suitable type of mattress platform according to the size of gaps between parts of mattress platform, which should not exceed 4 cm.
- ❖ If the mattress platform has not features to provide ventilation to the base of the mattress, then the mattress should be regularly lifted and the base allowed to dry/breath before further use. This will prevent moisture damage to the base cover and extend the life of the mattress.
- ❖ Continually loaded mattress has shorter lifetime period.

3.2 Contraindications

Mattresses are contraindicated for the following types of patients:

- patient with higher risk of pressure ulcers than the mattress is designed for
- patients with higher weight than mattress weight limit

3.3 Purpose of Use

Passive mattresses with increased antidecubitus effect. Their purpose is good weight distribution of lying person (patient) which contributes to prevention of pressure ulcers creation.

3.4 Conditions of Use

Mattresses must be used and stored in indoor environment where:

- The surrounding temperature is between +10 °C and +40 °C
- The relative humidity is between 30% and 75%
- The atmospheric pressure is between 795 hPa and 1060 hPa

4 Standards and Regulations

The mattresses comply with the following standards and directives:

- EN ISO 10993-5:2009
- EN ISO 10993-10:2010
- BS 7175:1989 (cover)
- BS 7177:2008+A1:2011 (mattress)
- EN 597-1:1994
- EN 597-2:1994

5 Variants

EffectaCare 10
 EffectaCare 20
 EffectaCare 20+ bariatric
 EffectaCare 20 paediatric
 PrimaCare 10
 PrimaCare 20
 CliniCare 10
 CliniCare 10 paediatric
 CliniCare 20
 CliniCare 20+ bariatric
 CliniCare 30

6 Delivery

⚠ CAUTION!

Damage to the mattress due to incorrect use!

- ➡ It is necessary to let mattress restore for at least 48 hours if the mattress is delivered in transport (compressed) packaging.
- ➡ Mattress cannot be in transport (compressed) packaging for more than 1 month. Mattress may lose its qualitative properties by violating this instruction.
- ➡ Mattress must be unpacked from the transport (compressed) packaging right after delivery.

The mattresses may be delivered in transport (compressed) packaging. This packaging is used for saving space when transporting and delivering mattress to the customer.





Shape of compressed packaging:




- Shape: cylinder
- Diameter of cylinder: 30-40 cm according to the height of mattress
- Package size: mattress width + approx. 5 cm

In case of mattress delivery in compressed packaging, restore the mattress as follows:

- ❖ Open the package cover.
 - ❖ Remove mattress from the package.
 - ❖ Place mattress on the flat surface.
 - ❖ Let the mattress restore to its original condition for approx. 48 hours.
- It is possible to start using mattress after 48 hours.

7 Description of Mattresses

Mattress name + picture	Mattress construction and description
<p>EffectaCare 10</p> 	 <p>The standard single-layer mattress designed for use in medical, nursing and home care. The mattress is designed for patients up to second degree risk* of developing pressure ulcers. Mattress is made of polyurethane foam with a washable, waterproof coating. The cover is sewn, made of biaxially elastic, breathable fabric that contains antibacterial ingredient. The zipper is located on long and short side.</p>
<p>EffectaCare 20</p> 	 <p>A single-layer mattress designed for use in medical, nursing and home care. The mattress is designed for patients up to second degree* risk of developing pressure ulcers. The mattress has two sided cutting. Cutting on the bottom of the mattress ensures that the mattress is very well adapted to the shape of the bed when positioning the patient. Cutting on top of the patient ensures patient comfort. Mattress is made of polyurethane foam with a washable, waterproof cover. The cover is sewn, made of biaxially elastic, breathable fabric that contains antibacterial ingredient. The zipper is located on long and short side.</p>
<p>PrimaCare 10</p> 	 <p>A preventive single-layer mattress designed for use in medical, nursing and home care. The mattress has longitudinal and transversal cutting on its whole surface. Cutting are concentrated in places that are most burdened by the weight of the patient. Given the technology described above, mattress is suitable for patients up to third degree* risk of developing pressure ulcers. Mattress is made of cold polyurethane foam with a washable, waterproof cover. The cover is sewn, made of four-way elastic (top cover and sides), breathable fabric that contains antibacterial ingredient. The zipper is located on long and short side.</p>
<p>PrimaCare 20</p> 	 <p>Preventive layered antidecubitus mattress suitable for use in health care, nursing and home care. Laying surface is made of cold polyurethane foam, bottom side is made of polyurethane foam. The mattress has two sided transversal cutting. Cutting on the bottom of the mattress ensures that the mattress is very well adapted to the shape of the bed when positioning the patient. Cutting on top of the mattress ensures patient comfort. Mattress is suitable for patients up to third degree* risk of developing pressure ulcers. The mattress is delivered with washable, waterproof cover. The cover is composed of all welded joints, made of four-way elastic (top cover and sides), breathable fabric that contains antibacterial ingredient. The zipper is located all around the periphery of mattress cover.</p>

Mattress name + picture	Mattress construction and description
<p data-bbox="161 259 304 286">CliniCare 10</p> 	 <p data-bbox="767 327 1426 831">Preventive three-layer antidecubitus mattress made of cold, polyurethane and memory foam with recommended use primarily in healthcare. Head and heel zones are made of perforated memory foam. The mattress is divided into three zones. The head zone, the heel zone and the zone for patient's body. Head zone and heel zone are composed of three foams different densities, hardness and profiling. In both zones are perforated cuts from memory foam. In the zone for the patient's body is transversal cutting. All of above mentioned unique technologies have been used to reach maximum patient comfort. Cutting on the bottom of the mattress ensures that the mattress is very well adapted to the shape of the bed when positioning the patient. Mattress edges are from stiffer foam, which provides greater stability for the patient and the necessary support when sitting on the edge of the bed. Given the technology described above, mattress is suitable for patients up to fourth degree* risk of developing pressure ulcers. The mattress is delivered with washable, waterproof cover. The cover is composed of all welded joints, made of four-way elastic (top cover and sides), breathable fabric that contains antibacterial ingredient. The zipper is located all around the periphery of mattress cover.</p>
<p data-bbox="161 848 304 875">CliniCare 20</p> 	 <p data-bbox="767 916 1426 1420">Preventive combined antidecubitus mattress made of cold and polyurethane foams with stiffed edges, laying surface made of visco elastic foam with high density weight and with seven zones of profilation. Mattress is recommended to use primarily in healthcare. The mattress is composed of three types of foams. The bottom foam layer and a stiffer foam layer, which is placed on the edge of mattress, are covered by a layer of seven zone profile from viscoelastic foam. All zones are deployed to ensure maximum support for the patient at all risky locations. Viscoelastic foam is perfectly adapted to each patient, which guarantees its superior comfort. Cutting on the bottom of the mattress ensures that the mattress is very well adapted to the shape of the bed when positioning the patient. Mattress edges are from stiffer foam, which provides greater stability for the patient and the necessary support when sitting on the edge of the bed. Given the technology described above, mattress is suitable for patients up to fourth degree* risk of developing pressure ulcers. The mattress is delivered with washable, waterproof cover. The cover is composed of all welded joints, made of biaxially elastic, breathable fabric that contains antibacterial ingredient. The zipper is located all around the periphery of mattress cover.</p>
<p data-bbox="161 1438 304 1464">CliniCare 30</p> 	 <p data-bbox="767 1505 1426 1951">Preventive combined antidecubitus mattress with THERMIC effect with recommended use primarily in healthcare. The mattress is composed of two types of foams. The lower layer of the mattress is made of classical foam with cutting on both sides of the foam. Due to cutting on the loading foam, mattress responds better when the patient lies down on the mattress. . Cutting on the bottom of the mattress ensures that the mattress is very well adapted to the shape of the bed when positioning the patient. The upper foam is viscoelastic with THERMIC effect, to ensure maximum comfort for the patient. THERMIC effect is the latest technological innovation that due to abilities of special foam, effectively directs the heat generated by the patient and thus contributes to significantly reduce the likelihood of developing pressure ulcers. Given the technology described above, mattress is suitable for patients up to fourth degree* risk of developing pressure ulcers. The mattress is delivered with washable, waterproof cover. The cover is composed of all welded joints, made of biaxially elastic, breathable fabric that contains antibacterial ingredient. The zipper is located all around the periphery of mattress cover.</p>

*with an appropriate turning and re-positioning regime

8 Using the Mattress

8.1 Turning and Rotating the Mattress

CAUTION!

Damage to the mattress due to incorrect mattress platform shapes and dimensions!

- Always consider shape and dimensions of the mattress platform for which the mattress will be used.

CAUTION!

Damage to the mattress due to not turning/rotating mattress in regular intervals!

- It is necessary to turn/rotate mattress in regular intervals to avoid permanent deformation of the mattress core and to extend service life of the mattress

WARNING!

Risk of injuring the patient due to incorrect placement of the patient on the mattress!

- Patient must be always laid on the surface with screen printing (e.g. numbers 1, 2, 3 & 4).

Mattresses may be single-sided or double-sided according to the mattress type. Single-sided mattresses have numbers 1 and 2 in the corners of top surface and can only be rotated not turned. Double-sided mattresses have numbers 1 and 2 in the corners of top surface and number 3 and 4 in corners of the bottom surface. Both surfaces of the mattress can be used as the mattress can be turned and rotated. When replacing the mattress cover on single sided types make sure that the protective flap on the cover zip is correctly positioned to prevent fluid entering via the zip as shown in Fig. 1.

To maintain performance of the mattress it is necessary to rotate the mattress or if the mattress is double-sided then also turn the mattress regularly. Recommended interval of rotating/turning is once a month. This is preventive measure against permanent deformation of the mattress core. When rotating or turning mattress, always grab mattress from bottom side (Fig. 2). Never manipulate with mattress by grabbing it by its seams. This may cause damaging the cover and risk of fluids penetrating into the mattress core.


Fig. 1 Correct placement of zip cover (cover over zip)

- Top cover
- Base cover


Fig. 2 Manipulation with mattress

8.1.1 Rotating

Rotating mattress (Fig. 3 – position 1) means exchanging the position of head and foot end of the mattress.

8.1.2 Turning

Turning mattress (Fig. 3 – position 2) means exchanging the top and bottom surface. Turning is only possible for double-sided mattresses. Mattress types can be found in chapter “Technical Specifications.”


Fig. 3 Rotating and turning of mattress

8.2 Mattress Distinctions

Single-sided mattress
(number 1 and 2 only)

Double-sided mattress

(top surface)

(bottom surface)


Fig. 4 Mattress types

ENGLISH

9 Technical Specifications

	EffectaCare 10	EffectaCare 20	PrimaCare 10
Material of Foam	polyurethane foam	polyurethane foam	polyurethane foam
Material of Cover Backing Coating	100% polyester 100% polyurethane	100% polyester 100% polyurethane	100% polyester 100% polyurethane
Hydrostatic head	>200 cm	>200 cm	>200 cm
Mattress shape	square	square	square
Cover	water resistant and vapour permeable antibacterial additive - silver		
Single / Double-sided	double-sided	single-sided	single-sided
Load limit	Load limit, weight and height of mattress differs according the model number of mattresses. For information about specific mattress please look at serial label located inside mattress cover.		
Weight of mattress			
Height of mattress			
Risk of pressure ulcers creation (depending on patient's condition)*	Up to risk II	Up to risk II	Up to risk III
Mattress type	standard, antidecubitus	preventive, antidecubitus	preventive, antidecubitus
Warranty of Foam*	1 year	1 year	2 years
Warranty of Cover**	1 year	1 year	2 years

	PrimaCare 20	CliniCare 10	CliniCare 20+30
Material of Foam	cold and polyurethane foam	cold, polyurethane and memory foam	viscoelastic, cold and polyurethane foam
Material of Cover Backing Coating	100% polyester 100% polyurethane	100% polyester 100% polyurethane	100% polyester 100% polyurethane
Hydrostatic head	>200 cm	>200 cm	>200 cm
Mattress shape	square	square	square
Cover	water resistant and vapour permeable antibacterial additive - silver		
Single / Double-sided	single-sided	single-sided	single-sided
Load limit	Load limit, weight and height of mattress differs according the model number of mattresses. For information about specific mattress please look at serial label located inside mattress cover.		
Weight of mattress			
Height of mattress			
Risk of pressure ulcers creation (depending on patient's condition)*	Up to risk III	Up to risk IV	Up to risk IV
Mattress type	preventive, antidecubitus	preventive, antidecubitus	preventive/combined antidecubitus
Warranty of Foam*	2 years	2 years	2 years
Warranty of Cover**	2 years	2 years	2 years

*with an appropriate turning and re-positioning regime

**while respecting the proper usage and cleaning instructions (according chapter "Cleaning/Disinfection")

General foam tolerances:

- Length: ± 10mm
- Width: ± 5mm
- Height: ± 5mm

10 Cleaning/Disinfection

⚠ CAUTION!
<p>Incorrect cleaning/disinfection can damage the mattress!</p> <ul style="list-style-type: none"> ➤ Do not use pressure or steam cleaners. ➤ Follow the instructions and observe the dosages recommended by the manufacturer. ➤ Ensure that disinfectants are selected and applied by qualified hygiene experts only.

10.1 General Guidance

For safe and gentle cleaning:

- ❖ Do not use any strong acids or alkalines, (optimum pH range 6 – 8. Do not exceed pH of 9).
- ❖ Only use detergents that are suitable for cleaning medical equipment.
- ❖ Do not use abrasive powders, steel wool, or other material and cleaning agents that might damage the mattress. Do not scrub mattress surface.
- ❖ Never use any corrosive or caustic detergents.
- ❖ Never use detergents that deposit calcium carbonate.
- ❖ Never use detergents with solvents that might affect the structure and consistency of the plastics (benzene, toluene, acetone etc.).
- ❖ Use only hospital-approved cleaners and observe local directives concerning infection control.
- ❖ Always rinse with water after cleaning and dry thoroughly before use.
- ❖ Observe local directives concerning infection control.

Mattress parts to be cleaned	Recommended Cleaning Agents (General cleaning)
Top Cover , Bottom Cover	Standard hospital detergents, Alcohol or Quaternary Ammonium based disinfectants, Chlorine based disinfectants containing up to 1000 ppm Chlorine, followed by rinsing with water and drying thoroughly before use.
	Decontamination: Blood spills/C-diff. etc
	Chlorine based disinfectants containing up to 10,000 ppm Chlorine. Dwell time on surface at 10,000 ppm of 5 minutes, followed by rinsing with water and drying thoroughly before use.
Mattress Core	Do not clean!

Due to the variety of laundry equipment, chemicals and conditions in use, customers should satisfy themselves through pre-testing. It is essential that cover be thoroughly rinsed and dried after all cleaning procedures and before storage or reuse. Wet or damp PU surfaces are more prone to mechanical damage than when dry.

As stated above, after application of a suitable cleaner, the surface must be rinsed with water and dried before use. (Even if the cleaner instructions say that this is not required). This prevents a build-up of chemicals on the mattress surface which could be reactivated during use and affect biocompatibility.

NOTE: Continued use of high concentration, chlorine-based disinfectants may significantly reduce the performance and the working life of a coated material.

Type of Cleaning	Parts to be cleaned
Routine Cleaning and Disinfection	■ external of mattress cover
Full Cleaning and Disinfection	■ external of mattress cover

10.2 Routine Cleaning and Disinfection

Cleaning the mattress:

- ❖ Check mattress cover top for any signs of damage or for liquid ingress.
- ❖ Replace or repair and completely disinfect mattress cover top if damaged. Also check if the mattress core is not contaminated. In case of core contamination, do not use the mattress and dispose the core.
- ❖ Leave mattress cover on mattress.
- ❖ Clean with 50 °C warm water with cleaning detergent.
- ❖ Rinse mattress with cold water.
- ❖ Let mattress air dry or wipe dry.
- ❖ Wipe mattress with disinfectant and rinse mattress with cold water.
- ❖ Let mattress dry or wipe dry.

10.3 Complete Cleaning and Disinfection

Cleaning Top/Bottom Cover:

Use standard hospital detergents, Alcohol based cleaners or Quaternary Ammonium based disinfectants. Suitable Chlorine based cleaners can be used at a concentration of 1000ppm. Stronger concentrations of chlorine can be used if required, (up to 10,000ppm), with a maximum dwell time of five minutes followed by rinsing with water and drying thoroughly before use..

After application of a suitable cleaner, the surface should be rinsed with water and dried before use. (Even if the cleaner instructions say that this is not required). This prevents a build up of chemicals on the mattress surface that could reactivate during use and affect biocompatibility.

Cleaning the mattress:

- ❖ Check mattress cover top and base for any signs of damage.
- ❖ Replace or repair and completely disinfect mattress cover top and base if damaged. Also check if the mattress core is not contaminated. In case of core contamination, do not use the mattress and dispose the core.
- ❖ Leave mattress cover on mattress.
- ❖ Clean all mattress cells and pipes with 50 °C warm water with cleaning detergent.
- ❖ Rinse mattress with cold water.
- ❖ Let mattress air dry or wipe dry.
- ❖ Wipe mattress with disinfectant.
- ❖ Rinse mattress with cold water.
- ❖ Let mattress air dry or wipe dry.

Machine washing of the top/base mattress covers:

- ❖ Remove cover (see Removing the Mattress Cover).
- ❖ If machine washing mattress top/base covers, the temperature should be raised during the wash cycle, to 65°C/149°F, for 10 -15 minutes, or 71°C/160° F, for 3 - 10 minutes, using hospital approved detergents and rinsing agents.
- ❖ Dry cover in tumble dryer at low temperature.

NOTE: *Maximum wash temperature 75°C/167°F.*

10.4 Mattress Core

The entire core of the mattress does not require any major cleaning. The core does not need disinfection. Once a month it is recommended to ventilate the mattress core (remove the mattress cover and leave the mattress core on ventilated area for 12 -24 hours). The mattress core cannot be washed by water or by disinfection.

10.5 Removing the Mattress Cover

- ❖ Carefully open zippers under side skirt of mattress cover.
- ❖ Remove top surface cover. Inspect cover and clean if necessary.
- ❖ Remove bottom surface cover (this only relates to mattresses with 360° zip).

After cleaning the mattress cover:

- ❖ Reinstall mattress cover by reversing the process described above.

NOTE: Cover of single-sided mattress must be put back on the foam core in the correct way (i.e. profiled side of the mattress must be uppermost, screen printed side of the cover must be uppermost) (see Description of Mattresses).

NOTE: Cover must be returned on the original core from which it was removed. Identification of mattresses cores can be found in chapter “Technical Specifications.”

10.6 Cover Cleaning and Disinfection Overview

Cover (Top & Base) Cleaning & Disinfection Requirements	EffectaCare 10	EffectaCare 20	PrimaCare 10
Machine washable cover	min. 71°C for 3 min.		
Tumble dry	low temperature 1 hour max		
1000 ppm chlorine	max exposure before rinsing 30 min.		
10,000 ppm chlorine	max exposure before rinsing 5 min.		

Cover (Top & Base) Cleaning & Disinfection Requirements	PrimaCare 20	CliniCare 10	CliniCare 20+30
Machine washable cover	min. 71°C for 3 min.		
Tumble dry	low temperature 1 hour max		
1000 ppm chlorine	max exposure before rinsing 30 min.		
10,000 ppm chlorine	max exposure before rinsing 5 min.		

11 Storage

When mattress is not in use:

- ❖ Clean and disinfect mattress.
- ❖ Let the mattress dry completely.
- ❖ Pack in suitable cover.
- ❖ Store in a place suitable for medical devices.

NOTE: Never store more than a 5 pieces of mattress on themselves for long time. Long time storage is considered storage longer than 14 days since the delivery date. Failure to comply this instruction may lead to the damage of the mattress.

12 Disposal

12.1 Environmental Protection

Linet® is aware of the important role that the protection of our environment plays for future generations.

The materials of this product are environmentally compatible. It does not contain hazardous substances on the basis of cadmium, mercury, asbestos, PCB or CFC.

All waste packaging left over from when the product is made operational is marked in accordance with the legal regulations on packaging. Sort the waste packaging, which is left over once the product is made operational by following the graphic symbols and deliver it to a person authorized to further utilize it.

The materials of the appliance are reusable. By reusing, material recycling or other forms of use of old appliances you give an important contribution to the protection of our environment.

12.2 Disposal

12.2.1 Within Europe

To dispose the appliance:

- ❖ When you dispose of your appliance do not put it into the household waste.
- ❖ Fines may be imposed under national regulations if an incorrect procedure for waste removal is used.

Protect your health and the environment. Thank you.

13 Warranty

Linet® will only be held responsible for the safety and reliability of products that are regularly serviced and used in accordance with the safety guidelines.

Should a serious defect arise that cannot be repaired during maintenance:

- ❖ Do not continue to use mattress.

Those products are covered by a 12-24 month warranty from the date of purchase. The length of warranty for each type of mattress and core is stated in chapter "Technical Specifications". The warranty covers all material and manufacturing-related failures and errors. Failures and errors caused by incorrect use and external effects are not covered. Justified complaints will be fixed free of charge during the warranty period. Proof of purchase, with the date of purchase, is required for all warranty service. Our standard terms and conditions apply.

14 EC Declaration of Conformity

EC CONFORMITY DECLARATION

Number PSEN0022
Version 01

Date and place of issue: 09. 07. 2015, Želečvice

Conformity declaration issued by:

Commercial name	Linnet spol. s r. o.
Registered address	Želečvice 5, 274 01 Slaný, Czech Republic
Reg. No.	00507814
Telephone	+420 312 576 111
Fax	+420 312 522 668

As the producer of the product - name (brand):	EffectaCare, CliniCare, PrimaCare		
Variants of the product:	EffectaCare 10	CliniCare 10	PrimaCare 10
	EffectaCare 20	CliniCare 20	PrimaCare 20
		CliniCare 30	
	4PP (Variants are specified in the technical documentation of the product).		
Description and function designation:	Passive preventive antidecubitus mattress for medical, nursing and home care. This EC conformity declaration also covers all applicable accessories.		
Classification of the product as the medical device:	Class I nonsterile, without measuring function, according to annex IX MDD 93/42/EEC – rule 1		

A) Declaration

I declare that the said product is safe under the conditions of common use in compliance with the instructions and that measures have been taken to ensure the conformity of all the products brought to market with basic requirements of directives related thereto, stated in paragraph B.

B) Fulfilled technical requirements of related regulations

This product's characteristics comply with the technical parameters related to it and stated in MDD 93/42/EEC which stipulates the technical parameters for healthcare products.

C) Means of assessing conformity

Conformity was assessed by the procedure stated MDD 93/42/EEC, Annex VII

D) Used standards for product conformity assessment

The said product fulfills the requirements of these harmonized technical standards: ISO 10993-5:2009, ISO 10993-10:2010 and non-harmonized technical standards: BS 7175:1989. which were used for assessing of conformity

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Copyright © LINET, s.r.o., 2017